



Test a set of external controls with each new shipment received and once for each untrained operator. Further controls may be tested in order to conform with local, state and/or federal regulations, accrediting groups, or your lab's standard Quality Control procedures. An operator should confirm a valid test prior to interpreting results. A valid test is a green check mark present as well as the presence of a results valid purple spot. For more information, please review Instructions for Use and Quick Reference Guide.

Note: If the positive or negative external controls fail, repeat with a new Visby Medical Sexual Health Test. If the repeat test fails, please contact: Visby Medical Customer Support at 1-833-468-4729 (1-833-GoVisby).

	Device Information		Control Information		Internal Controls				External Controls		
Date	Lot #	Exp Date	Lot #	Exp Date	#1 Valid Electronic Control (Green Check Mark-Present or Absent)		#2- Results Valid Spot (Spot Development-Present or Absent)		#3-5 Valid Positive Control? (Spot development for all pathogens present)	#3-5 Valid Negative Control (No spot development for pathogens present)	Tech Initials
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	

Date	Device Information		Control Information		Internal Controls				External Controls		Tech Initials
	Lot #	Exp Date	Lot #	Exp Date	#1 Valid Electronic Control (Green Check Mark-Present or Absent)		#2- Results Valid Spot (Spot Development-Present or Absent)		#3-5 Valid Positive Control? (Spot development for all pathogens present)	#3-5- Valid Negative Control (No spot development for pathogens present)	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	
/ /		/ /		/ /	Present	Absent	Present	Absent	Present	Absent	

Comments